

POCT

Immunofluorescence Analyzer



Immunofluorescent Analyzer is an in-vitro diagnostic device for the quantitative analysis of test kits detecting a wide range of antigens, antibodies, hormones, and other proteins.



Advantages

- Maximized sensitivity
- Self-Calibration
- Android system
- LIS/HIS connectivity
- Small size, easy to carry

PERFORMANCE PARAMETER

Product Name	Immunofluorescence Analyzer
Model	VMFIA1001
Size	245*270*160 mm
Display	8 inches
Methodology	Time-resolved fluorescence Immunoassay
Weight	<2 kg
Sample type	Whole Blood/Serum/Plasma/Urine/etc.
Test channel	Single channel, automatic
Test items	Up to 63
Test mode	Single cassette multi-project (≤4 items)
Test speed	Single cassette < 10s
Internet	Serial Port, USB, WLAN
Print	Built-in thermal printer
System	Android
Storage	> 10000 test records
Language	English, French, Spanish, Portuguese, etc.
Card code	QR code

ORDER INFORMATION

Category	Cat. No	Name	Specimen	Time	Kit Size	CE Status
Hormone	VMPO15	β-HCG Rapid Test Kit	S/P/WB	15 mins.	25 T	CE
	VMPO36	Progesterone (PROG) Rapid Test Kit	S/P	15 mins.	25 T	CE
	VMPO37	Prolactin (PRL) Rapid Test Kit	S/P/WB	15 mins.	25 T	CE
	VMPO38	Luteinizing Hormone (LH) Rapid Test Kit	S/P/WB	15 mins.	25 T	CE
	VMPO39	Follicle Stimulating Hormone (FSH) Rapid Test Kit	S/P/WB	15 mins.	25 T	CE
	VMPO41	Testosterone (TES) Rapid Test Kit	S/P	15 mins.	25 T	CE
	VMPO42	Anti-Mullerian Hormone (AMH) Rapid Test Kit	S/P/WB	15 mins.	25 T	CE
	VMPO43	Estradiol (E2) Rapid Test Kit	S/P/WB	15 mins.	25 T	CE
COVID-19	VSCD03	SARS-CoV-2 Rapid Antigen Test (FIA)	NP	10 mins.	25 T	CE
	VSCD04	SARS-CoV-2 NAb Test (FIA)	S/P/WB	15 mins.	25 T	CE
Inflammation	VMPO29	Serum Amyloid A Determination Kit (SAA Determination Kit)	S/P/WB	10 mins.	25 T	CE
	VMPO09	Procalcitonin (PCT) Rapid Test Kit	S/P/WB	15 mins.	25 T	CE
	VMPO11	C-reactive Protein (hsCRP&CRP) Rapid Test Kit	S/P/WB	3 mins.	25 T	CE
	VMPO34	Interleukin 6 (IL-6) Rapid Test Kit	S/P/WB	15 mins.	25 T	CE
Thyroid	VMPO01	Thyroid-Stimulating Hormone (TSH) Rapid Test Kit	S/P/WB	15 mins.	25 T	CE
	VMPO02	Triiodothyronine (T3) Rapid Test Kit	S/P/WB	15 mins.	25 T	CE
	VMPO03	Thyroxine (T4) Rapid Test Kit	S/P/WB	15 mins.	25 T	CE
	VMPO13	Free T3 (FT3) Rapid Test Kit	S/P/WB	15 mins.	25 T	CE
	VMPO14	Free T4 (FT4) Rapid Test Kit	S/P/WB	15 mins.	25 T	CE
Tumor Markers	VMPO07	Prostate Specific Antigen (PSA) Rapid Test Kit	S/P/WB	15 mins.	25 T	CE
Bone	VMPO31	25-hydroxy Vitamin D (25-OHVD) Rapid Test Kit	S	10 mins.	25 T	CE
Cardiac Marker	VMPO04	Cardiac Troponin I (cTnI) Rapid Test Kit	S/P/WB	15 mins.	25 T	CE
	VMPO05	N-terminal Prohormone of Brain Natriuretic Peptide (NT-pro BNP) Rapid Test Kit	S/P/WB	15 mins.	25 T	CE
	VMPO08	D-Dimer (DD) Rapid Test Kit	S/P/WB	5 mins.	25 T	CE
	VMPO12	Myoglobin (Myo) Rapid Test Kit	S/P/WB	15 mins.	25 T	CE
	VMPO45	BNP Rapid Test Kit	P/WB	15 mins.	25 T	CE
	VMPO46	Cardiac Troponin T (cTnT) Rapid Test Kit	S/P/WB	15 mins.	25 T	CE
	VMPO18	Creatine Kinase Isoenzymes MB (CK-MB) Rapid Test Kit	S/P/WB	15 mins.	25 T	CE
Diabetes	VMPO06	Hemoglobin A1c (HbA1c) Rapid Test Kit	WB	5 mins.	25 T	CE