

VNFOR

Viral Nucleic Acid Transport Medium

VNFOR Viral Nucleic Acid Transport Medium is an extraction buffer used to extract viral nucleic acids from clinical samples or media samples to obtain viral nucleic acids for use in nucleic acid amplification tests.

Advantages

- Contamination-preventing content
- Flexibility in transport temperature
- Suitable for further testing such as RT-PCR
- Crack-resistant, leak-proof, and stand-up tubes
- Fast, simple, and inexpensive



Ref. No	Name	Storage	Kit Size	CE Status
VNF01	Viral Nucleic Acid Transport Medium	2°C~ 25°C	100 T	CE