

RapidFor™



EASY



ACCURATE



RAPID




RESULT IN 15 MINUTES

SARS-CoV-2 Antigen Antigen Test


The RapidFor™ Saliva Rapid Antigen Test Kit uses colloidal label as a tool for SARS-CoV-2 antigen detection in the collected saliva, saliva swab specimen.

Avantages

- Easy to use
- Pain-free, no swab
- No extra lab equipment
- Variant Detection [B.1.1.7 (Alpha), B.1.351 (Beta), E484K mutation, B1.617.2 (Delta), P.1 (Gamma), C.37 (Lambda), B.1.1.529 (Omicron)]
- Early phase detection
- Suitable for kids




How to Use



Box Content

25x Test Cassette
25x Extraction Buffer
25x Collecting Apparatus
1x Instruction For Use

Interpretation



Cat. No	Target	Specimen	Kit Size	Carton Size	CE Durumu
VSCD05	SARS-CoV-2 Antigen	Saliva	25T	58x45x43cm	CE