

FOB TRANSFERRIN RAPID TEST KIT




FOB + Transferrin Rapid Test Kit

FOB + Transferrin Rapid Test Kit is designed for the in vitro qualitative detection and differentiation of fecal occult blood (FOB) and transferrin in human fecal samples.

Advantages

- Assists the diagnosis of gastrointestinal bleeding disorders
- Results are available within 15 minutes
- Easy to use procedure
- Visually interpretation
- No need for extra lab analyzer


Test Procedure



Box Content

- 25x Test Cassette
- 25x Extraction Buffer
- 1x Instruction For Use

Interpretation



Cat. No	Target	Specimen	Kit Size	CE Status
VSCD70	Human Occult Blod	F	25T	CE